


# Webly - Lite Ankle Orthosis


**Warnings and Instructions: Review carefully, proper application is required**

**⚠ Warning:** This device will not prevent or eliminate risk of injury. Do Not Over-tighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.


**Care:** Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. The Webly-Lite Orthosis is designed to be worn over an athletic sock. Position your heel into the hole at the rear and bottom of brace. Tighten the laces snugly and tie.



2. Bring the inside strap forward, across the top of the foot. Wrap it under the foot, continuing under the heel, and pull the strap upward securing the hook fastener onto the loop material on the inside of the ankle.



3. Bring the outside strap forward, across the top of the foot. Wrap it under the foot, continuing under the heel, and pull the strap upward securing the hook fastener onto the loop material on the outside of the ankle.



4. Stand, putting full weight on your ankle, pull up on both straps simultaneously to gain proper alignment and tightness. Attach hook closures.

**HELY & WEBER**

© 2007 Weber Orthopedic Inc.  
DBA Hely & Weber All rights reserved.  
1185 East Main St., Santa Paula, California 93060  
California: 800-221-5465, National: 800-654-3241  
International: 805-525-4244, fax: 805-933-2348  
U.S. fax: 800-559-5975, [www.hely-weber.com](http://www.hely-weber.com)

**CE**

**EC REP**

MDSS GmbH  
Schiffgraben 41  
30175 Hannover, Germany

Rev A 11/08