

455



Removable, malleable palmar stay



Removable, malleable dorsal stay

TITAN™ THUMB ORTHOSIS

- Constructed of soft, comfortable orthopedic felt laminated to brush nylon exterior
- The adjustable lacing pull easily fits a wide range of wrist sizes
- Oblique thumb stay is malleable to allow for desired thumb position and rigid enough for excellent thumb stability
- One size right and one size left significantly reduces inventory
- Provides abduction control, angle adjustability, and easy reapplication


Order No. Description / Options
455 Titan™ Thumb Orthosis
459 Long Titan™ Thumb Orthosis
Specify Right or Left




459



TITAN™ THUMB SIZING:
 Wrist circumference 4.5" to 9.5"



The oblique thumb stay can be bent to allow desired thumb positioning



On small hands, the eyelets will clamshell over the laces for a smooth overlap



The brace remains below the palmar crease and allows full finger function while securing the thumb

