

450



The brace fits below the palmar crease and allows full finger function

Anatomically Contoured Aluminum Palmar Stay



TITAN™ WRIST - LACING ORTHOSIS

- Constructed of soft, comfortable orthopedic felt laminated to brush nylon exterior
- Universal size reduces stocking needs
- Adjustable pull and lacing system handles a wide range of wrist sizes
- Anatomically contoured palmar stay is malleable for customization and removable for "step-down" rehabilitation

Order No. Description / Options
450 Titan™ Wrist
452 Titan™ Wrist & Forearm
Specify Right or Left


TITAN™ WRIST SIZING:

Wrist Circumference 4.5" to 9.5"



452

Dorsal stay pod can be moved proximal or distal for desired level of dorsiflexion restriction



Dorsal range of motion can be controlled by changing the position of the dorsal stay pod. Move the stay pod distal to restrict dorsiflexion. Move the stay pod proximal for less dorsiflexion restriction.



On small hands, the eyelets will clamshell over the laces for a smooth overlap



One pull design makes application simple