


# Cervical Collar


619  
629  
639

## Warnings and Instructions: Review carefully, proper application is required

**⚠ Warning:** This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.  
Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.



1. Place under neck, centering the collar under the chin and leaving the head in a neutral position. **DO NOT HYPEREXTEND OR BEND THE NECK BACKWARDS.** Wrap collar around the neck.
2. Apply the contact closure to secure the collar.



3. Finished Application.

**HELY&WEBER**

© 2014 Weber Orthopedics Inc.  
DBA Hely & Weber all rights reserved

1185 E Main St., Santa Paula, California 93060  
California: 800-221-5465, National: 800-654-3241  
International: 805-525-4244, fax: 805-933-2348  
U.S. fax: 800-559-5975, www.hely-weber.com

**CE**  
[EC] [REP]

MDSS GmbH, Schiffgraben 41  
30175 Hannover, Germany


Rev B 2/14

# Cervical Collar


619  
629  
639

## Warnings and Instructions: Review carefully, proper application is required

**⚠ Warning:** This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.  
Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.



1. Place under neck, centering the collar under the chin and leaving the head in a neutral position. **DO NOT HYPEREXTEND OR BEND THE NECK BACKWARDS.** Wrap collar around the neck.
2. Apply the contact closure to secure the collar.



3. Finished Application.

**HELY&WEBER**

© 2014 Weber Orthopedics Inc.  
DBA Hely & Weber all rights reserved

1185 E Main St., Santa Paula, California 93060  
California: 800-221-5465, National: 800-654-3241  
International: 805-525-4244, fax: 805-933-2348  
U.S. fax: 800-559-5975, www.hely-weber.com

**CE**  
[EC] [REP]

MDSS GmbH, Schiffgraben 41  
30175 Hannover, Germany

Rev B 2/14